



Beyond limits

EverGo portable oxygen concentrator offers up to 8 hours of battery life, infinite possibilities.

PHILIPS

RESPIRONICS



Endless



EverGo — the right choice for your business and your patients

Discover the POC that makes business sense for you

- Provides a new retail opportunity
- Supports patient travel programs
- Features are attractive to patients
- Differentiates your oxygen business to physicians
- Helps control rising costs by pairing with EverFlo for a complete oxygen generating system

Eight hours of battery life* take your patients beyond limits

- 8 hours of battery life is a full day of freedom for your patients
- 24/7 global support network lets you rest assured when your patients travel outside your area
- Federal Aviation Administration (FAA)-approved for all flights departing and landing in the United States**

possibilities

*Based on 20 breaths per minute and a device setting of two.

** Airline policies for traveling with an oxygen concentrator vary. Patients need to check with their airlines before flying with EverGo.

EverGo all



inclusive

The EverGo Portable Oxygen Concentrator ships with everything included. There is no need for providers or oxygen users to purchase additional accessories

The Freedom Series is a family of products designed to deliver more freedom – freedom from worry and freedom from complicated maintenance. The Freedom Series frees you and your patients from some of the more demanding aspects of oxygen therapy, so you both can go on with your lives.



- Integrated EverGo carrying case
- 2 rechargeable lithium ion batteries that provide up to 8 hours of battery life*
- AC and DC power cords allow for charging at home or on the go
- Fold up, lightweight wheeled cart
- Convenient accessory case
- Lightweight (<10 lbs) non-medical design
- Easy-to-use touch screen display controls all functions and clearly shows battery life and setting
- Quiet at 43 dBA
- Low maintenance – one easy-to-remove patient filter
- Easy access to batteries makes them easy to change and charge – plus EverGo can be simultaneously charged and operated
- Large oxygen capacity (1050 ml/min)


*Other portable systems may require the purchase of as many as three extra batteries to match the operating time of EverGo's two included batteries.

Six dosage settings let you customize therapy to your patients' needs

Setting	Dosage
1	12 ml
2	23 ml
3	35 ml
4	47 ml
5*	59 ml*
6*	70 ml*

*Based on a respiratory rate of 15 breaths/minute

Portable oxygen concentrator battery life comparison



Product specifications

AC power	Universal input: 100-240 VAC, 50-60 Hz, 1.0 Amp Output: 18VDC, 5 Amp
DC power	12-18VDC, 5 Amp
Power consumption	110 W at 100VAC (charging) 120 W at 240 VAC (charging); 70 W (not charging)
Oxygen concentration	90% +/- 3%
Oxygen delivery pulse dose	Settings 1-6 in 0.5 increments
Weight	10 lbs (4.5 kg) with 2 batteries, 8.5 lbs (3.9 kg) with 1 battery
Dimensions	Length: 12 in (30.5 cm) Width: 6 in (15.2 cm) Height: 8.5 in (21.6 cm)
Battery duration	2 batteries = 8 hours (setting of 2, 20 bpm) 1 battery = 4 hours (setting of 2, 20 bpm)
Battery recharge time	2 hours – 1 battery 4 hours – 2 batteries (based on fully discharged battery)
Sound level	43 dBA typical (setting of 2, 20 bpm)
Filters	Gross particle, inlet filter, micro-felt filter (bacteria)
Cannula length	30 ft maximum
Operating temperature	41°F to 104°F (5°C to 40°C)
Storage/transport humidity	-4°F to 140°F (-20°C to 60°C) up to 95% RH non-condensing
Operating humidity	Up to 95% relative humidity
Outlet pressure	3 PSI
Operating altitude	Up to 8000 ft
Oxygen purity alarms	Low oxygen at 83%
Alarms and indicators	Low battery at 17 minutes life left; battery status on touch screen and on battery; no breath alarm after 2 minutes; high breath rate alarm
Oxygen capacity	1050 ml per minute (pulse dosed)
Pulse dose settings	1 - 12 ml, 2 - 23 ml, 3 - 35 ml, 4 - 47 ml,* 5 - 59 ml,* 6 - 70 ml

*Based on a respiratory rate of 15 breaths/minute

EverGo accessories

Ordering information

Basic system includes: EverGo portable oxygen concentrator, carrying case (choice of color), two rechargeable lithium ion batteries, mobile cart, DC and AC power cords, accessory case and user manual.

Concentrator systems	Part number
EverGo system with black carrying case	900-000
EverGo system with navy blue carrying case	900-002
EverGo system with black carrying case (Canada)	900-010
EverGo system with navy blue carrying case (Canada)	900-012
Accessories	
Air inlet filter (6 pack)	900-150
EverGo black carrying case	900-200
EverGo blue carrying case	900-202
Rechargeable li ion battery	900-102
Accessory case	900-250
Battery charger/recalibrator	900-103
Cart	900-100
DC power cord	900-101
AC power supply	900-105
AC power cord (North America)	900-104
User manual	900-098
Reusable ship case	900-205
Extended warranty	
EverGo 1-year extended warranty	1038653
EverGo 2-year accident protection	1038570
EverGo 3-year accident protection (includes 1-year extended warranty)	1038648
Maintenance plans	
Deluxe maintenance plan	1038651
Premium maintenance plan	1038652



Black carrying case
900-200



Rechargeable li ion battery
900-102



Accessory case
900-250



Battery charger/recalibrator
900-103



Cart
900-100



DC power cord
900-101



User manual
900-098

Philips Healthcare is part of
Royal Philips Electronics

How to reach us

www.philips.com/healthcare
healthcare@philips.com
fax: +31 40 27 64 887

Asia
+852 2821 5888

Europe, Middle East, Africa
+49 7031 463 2254

Latin America
+55 11 2125 0744

North America
+1 425 487 7000
800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)

Respironics Asia Pacific
+65 6298 1088

Respironics Australia
+61 (2) 9666 4444

Respironics Deutschland
+49 8152 93 06 0

Respironics Europe, Middle East, Africa
+33 1 47 52 30 00

Respironics France
+33 2 51 89 36 00

Respironics Italy
+39 03 62 63 43 1

Respironics Sweden
+46 8 120 45 900

Respironics Switzerland
+41 6 27 45 17 50

Respironics United Kingdom
+44 800 1300 845

www.philips.com/respironics

Please visit <http://evergo.respironics.com>

Respironics and EverGo are trademarks of Respironics, Inc. and its affiliates.



© 2009 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.
CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Geyer WMB 8/11/09 MCI 4102280 PN 1042845